



PSD3 rabbit pAb

Cat No.:ES18988

For research use only

Overview

Product Name	PSD3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PSD3 AA range: 116-166
Specificity	This antibody detects endogenous levels of human PSD3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	PH and SEC7 domain-containing protein 3 (Epididymis tissue protein Li 20mP) (Exchange factor for ADP-ribosylation factor guanine nucleotide factor 6) (Hepatocellular carcinoma-associated antigen 67) (
Gene Name	PSD3 EFA6R HCA67 KIAA0942
Cellular localization	Cell membrane . Cell projection, ruffle membrane . Cell junction, synapse, postsynaptic density . In interphase associated with the plasma membrane, in particular with membrane ruffling regions. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23362
Human Swiss-Prot Number	Q9NYI0
Alternative Names	
Background	

