

CRFM7 rabbit pAb

Cat No.:ES18989

For research use only

Overview

Product Name	CRFM7 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CRFM7 AA range: 124-174
Specificity	This antibody detects endogenous levels of human CRFM7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	CHRNA7-FAM7A fusion protein (CHRNA7-DR1) (D-10)
Gene Name	CHRFAM7A
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	89832
Human Swiss-Prot Number	Q494W8
Alternative Names	
Background	

