



ST1A3 rabbit pAb

Cat No.:ES18990

For research use only

Overview

Product Name	ST1A3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ST1A3 AA range: 227-277
Specificity	This antibody detects endogenous levels of human ST1A3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Sulfotransferase 1A3/1A4 (ST1A3/ST1A4) (EC 2.8.2.1) (Aryl sulfotransferase 1A3/1A4) (Catecholamine-sulfating phenol sulfotransferase) (HAST3) (M-PST) (Monoamine-sulfating phenol sulfotransferase) (Pla
Gene Name	SULT1A3 STM; SULT1A4
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	445329
Human Swiss-Prot Number	P50224
Alternative Names	
Background	

