



CGNL1 rabbit pAb

Cat No.:ES18994

For research use only

Overview

Product Name	CGNL1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CGNL1 AA range: 394-444
Specificity	This antibody detects endogenous levels of human CGNL1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Cingulin-like protein 1 (Junction-associated coiled-coil protein) (Paracingulin)
Gene Name	CGNL1 JACOP KIAA1749
Cellular localization	Cell junction, tight junction . Localizes to the apical junction complex composed of tight and adherens junctions (PubMed:22891260). In the liver and kidney, it is also found along non-junctional actin filament bundles in addition to the apical junction (By similarity). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84952
Human Swiss-Prot Number	Q0VF96
Alternative Names	
Background	

