



CX04A rabbit pAb

Cat No.:ES18995

For research use only

Overview

Product Name	CX04A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CX04A AA range: 15-65
Specificity	This antibody detects endogenous levels of human CX04A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Protein CXorf40A (Endothelial-overexpressed lipopolysaccharide-associated factor 1)
Gene Name	CXorf40A CXorf40 EOLA1
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	91966
Human Swiss-Prot Number	Q8TE69
Alternative Names	
Background	

