

## HS3S1 rabbit pAb

Cat No.: ES18996

For research use only

## Overview

Product Name HS3S1 rabbit pAb

Host species Rabbit Applications WB

Species Cross-Reactivity Human; Mouse; Rat Recommended dilutions WB 1:1000-2000

Immunogen Synthesized peptide derived from human HS3S1 AA

range: 201-251

**Specificity** This antibody detects endogenous levels of human

**HS3S1** 

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20℃. Avoid repeated freeze-thaw cycles.

Protein Name Heparan sulfate glucosamine 3-O-sulfotransferase 1

(EC 2.8.2.23) (Heparan sulfate D-glucosaminyl 3-O-sulfotransferase 1) (3-OST-1) (Heparan sulfate

3-O-sulfotransferase 1) (h3-OST-1)

Gene Name HS3ST1 3OST 3OST1
Cellular localization Golgi apparatus lumen .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 9957 Human Swiss-Prot Number 014792

**Alternative Names** 

+86-27-59760950

**Background** 

