

THUM3 rabbit pAb

Cat No.:ES18997

For research use only

Overview

Product Name	THUM3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human THUM3 AA range: 410-460
Specificity	This antibody detects endogenous levels of human THUM3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	THUMP domain-containing protein 3
Gene Name	THUMPD3
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	25917
Human Swiss-Prot Number	Q9BV44
Alternative Names	
Background	

