



F177A rabbit pAb

Cat No.:ES19000

For research use only

Overview

Product Name	F177A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human F177A AA range: 112-162
Specificity	This antibody detects endogenous levels of human F177A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Protein FAM177A1
Gene Name	FAM177A1 C14orf24
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	283635
Human Swiss-Prot Number	Q8N128
Alternative Names	
Background	

