

# CCD34 rabbit pAb

Cat No.:ES19003

For research use only

## Overview

Product Name	CCD34 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CCD34 AA range: 249-299
Specificity	This antibody detects endogenous levels of human CCD34
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Coiled-coil domain-containing protein 34 (Renal carcinoma antigen NY-REN-41)
Gene Name	CCDC34 RAMA3 L15
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	91057
Human Swiss-Prot Number	Q96HJ3
Alternative Names	
Background	

