

# WDR65 rabbit pAb

Cat No.:ES19004

For research use only

## Overview

Product Name	WDR65 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human WDR65 AA range: 457-507
Specificity	This antibody detects endogenous levels of human WDR65
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	WD repeat-containing protein 65
Gene Name	WDR65
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	149465
Human Swiss-Prot Number	Q96MR6
Alternative Names	
Background	

