



MYCT1 rabbit pAb

Cat No.:ES19006

For research use only

Overview

Product Name	MYCT1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MYCT1 AA range: 114-164
Specificity	This antibody detects endogenous levels of human MYCT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Myc target protein 1 (Myc target in myeloid cells protein 1)
Gene Name	MYCT1 MTLC MTMC1
Cellular localization	Nucleus . Expressed in nuclei of hepatocellular carcinoma cell line BEL-7402 cells.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	80177
Human Swiss-Prot Number	Q8N699
Alternative Names	
Background	

