

SPAT4 rabbit pAb

Cat No.:ES19007

For research use only

Overview

Product Name	SPAT4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SPAT4 AA range: 62-112
Specificity	This antibody detects endogenous levels of human SPAT4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Spermatogenesis-associated protein 4 (Testis and spermatogenesis cell-related protein 2) (Testis spermatocyte apoptosis-related gene 2 protein)
Gene Name	SPATA4 TSARG2
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	132851
Human Swiss-Prot Number	Q8NEY3
Alternative Names	
Background	

