

AGAP5 rabbit pAb

Cat No.:ES19010

For research use only

Overview

Product Name	AGAP5 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human AGAP5 AA range: 417-467
Specificity	This antibody detects endogenous levels of human AGAP5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Arf-GAP with GTPase, ANK repeat and PH domain-containing protein 5 (AGAP-5) (Centaurin-gamma-like family member 2)
Gene Name	AGAP5 CTGLF2
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	729092
Human Swiss-Prot Number	A6NIR3
Alternative Names	
Background	

