



FCGRC rabbit pAb

Cat No.:ES19011

For research use only

Overview

Product Name	FCGRC rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FCGRC AA range: 200-250
Specificity	This antibody detects endogenous levels of human FCGRC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Putative high affinity immunoglobulin gamma Fc receptor IC (IgG Fc receptor IC) (Fc-gamma RIC) (FcRIC) (hFcgammaRIC)
Gene Name	FCGR1C IGFRC
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	A6NKC4
Alternative Names	
Background	

