



NECP1 rabbit pAb

Cat No.:ES19013

For research use only

Overview

Product Name	NECP1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NECP1 AA range: 109-159
Specificity	This antibody detects endogenous levels of human NECP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Adaptin ear-binding coat-associated protein 1 (NECAP endocytosis-associated protein 1) (NECAP-1)
Gene Name	NECAP1
Cellular localization	Cytoplasmic vesicle, clathrin-coated vesicle membrane . Cell membrane . Colocalizes with AP-2 at the plasma membrane. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	25977
Human Swiss-Prot Number	Q8NC96
Alternative Names	
Background	

