



COMD5 rabbit pAb

Cat No.:ES19014

For research use only

Overview

Product Name	COMD5 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human COMD5 AA range: 47-97
Specificity	This antibody detects endogenous levels of human COMD5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	COMM domain-containing protein 5 (Hypertension-related calcium-regulated gene protein) (HCaRG)
Gene Name	COMMD5 HT002
Cellular localization	Cytoplasm . Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	28991
Human Swiss-Prot Number	Q9GZQ3
Alternative Names	
Background	

