



HMN3 rabbit pAb

Cat No.:ES19022

For research use only

Overview

Product Name	HMN3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HMN3 AA range: 1-24
Specificity	This antibody detects endogenous levels of human HMN3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Humanin-like protein 3 (HN3) (MT-RNR2-like protein 3)
Gene Name	MTRNR2L3
Cellular localization	Secreted . Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	100462983
Human Swiss-Prot Number	POCJ70
Alternative Names	
Background	

