



KLF9 rabbit pAb

Cat No.:ES19024

For research use only

Overview

Product Name	KLF9 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KLF9 AA range: 49-99
Specificity	This antibody detects endogenous levels of human KLF9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Krueppel-like factor 9 (Basic transcription element-binding protein 1) (BTE-binding protein 1) (GC-box-binding protein 1) (Transcription factor BTEB1)
Gene Name	KLF9 BTEB BTEB1
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	687
Human Swiss-Prot Number	Q13886
Alternative Names	
Background	

