



SPAT1 rabbit pAb

Cat No.:ES19027

For research use only

Overview

Product Name	SPAT1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SPAT1 AA range: 369-419
Specificity	This antibody detects endogenous levels of human SPAT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Spermatogenesis-associated protein 1 (Sperm-specific protein SP-2)
Gene Name	SPATA1
Cellular localization	Cytoplasmic vesicle, secretory vesicle, acrosome .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q5VX52
Alternative Names	
Background	

