

RET1G rabbit pAb

Cat No.: ES19031

For research use only

Overview

Product Name RET1G rabbit pAb

Host species Rabbit Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1:1000-2000

Immunogen Synthesized peptide derived from human RET1G AA

range: 192-242

Specificity This antibody detects endogenous levels of human

RET1G

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Retinoic acid early transcript 1G protein

Gene Name RAET1G

Cellular localization [Isoform 1]: Cell membrane; Lipid-anchor,

GPI-anchor . Endoplasmic reticulum . Mainly found

intracellularly. .; [Isoform 2]: Secreted .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/ml

Observed band

Human Gene ID 353091 Human Swiss-Prot Number Q6H3X3

Alternative Names

+86-27-59760950

Background

