



CBLC rabbit pAb

Cat No.:ES19036

For research use only

Overview

Product Name	CBLC rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CBLC AA range: 219-269
Specificity	This antibody detects endogenous levels of human CBLC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Signal transduction protein CBL-C (RING finger protein 57) (SH3-binding protein CBL-3) (SH3-binding protein CBL-C)
Gene Name	CBLC CBL3 RNF57
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23624
Human Swiss-Prot Number	Q9ULV8
Alternative Names	
Background	

