

TEX35 rabbit pAb

Cat No.:ES19037

For research use only

Overview

Product Name	TEX35 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TEX35 AA range: 96-146
Specificity	This antibody detects endogenous levels of human TEX35
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Testis-expressed sequence 35 protein
Gene Name	Tex35 C1orf49
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84066
Human Swiss-Prot Number	Q5T0J7
Alternative Names	
Background	

