



EX3L2 rabbit pAb

Cat No.:ES19038

For research use only

Overview

Product Name	EX3L2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human EX3L2 AA range: 271-321
Specificity	This antibody detects endogenous levels of human EX3L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Exocyst complex component 3-like protein 2 (HBV X-transactivated gene 7 protein) (HBV XAg-transactivated protein 7)
Gene Name	EXOC3L2 XTP7
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	90332
Human Swiss-Prot Number	Q2M3D2
Alternative Names	
Background	

