



DC8L2 rabbit pAb

Cat No.:ES19041

For research use only

Overview

Product Name	DC8L2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DC8L2 AA range: 458-508
Specificity	This antibody detects endogenous levels of human DC8L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	DDB1- and CUL4-associated factor 8-like protein 2 (WD repeat-containing protein 42C)
Gene Name	DCAF8L2 WDR42C
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	POC7V8
Alternative Names	
Background	

