



SRSF2 rabbit pAb

Cat No.:ES19042

For research use only

Overview

Product Name	SRSF2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SRSF2 AA range: 15-65
Specificity	This antibody detects endogenous levels of human SRSF2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Serine/arginine-rich splicing factor 2 (Protein PR264) (Splicing component, 35 kDa) (Splicing factor SC35) (SC-35) (Splicing factor, arginine/serine-rich 2)
Gene Name	SRSF2 SFRS2
Cellular localization	Nucleus . Nucleus, nucleoplasm . Nucleus speckle . Phosphorylation by SRPK2 provokes its redistribution from the nuclear speckle to nucleoplasm.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	6427
Human Swiss-Prot Number	Q01130
Alternative Names	
Background	

