



LUZP4 rabbit pAb

Cat No.:ES19044

For research use only

Overview

Product Name	LUZP4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LUZP4 AA range: 42-92
Specificity	This antibody detects endogenous levels of human LUZP4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Leucine zipper protein 4 (Cancer/testis antigen 28) (CT-28) (CT28) (Tumor antigen HOM-TES-85)
Gene Name	LUZP4
Cellular localization	Nucleus . Cytoplasm . In nuclear speckles. Relocalizes to the cytoplasm during cell division. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	51213
Human Swiss-Prot Number	Q9P127
Alternative Names	
Background	

