



ASCL4 rabbit pAb

Cat No.:ES19045

For research use only

Overview

Product Name	ASCL4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ASCL4 AA range: 81-131
Specificity	This antibody detects endogenous levels of human ASCL4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Achaete-scute homolog 4 (ASH-4) (hASH4) (Achaete-scute-like protein 4) (Class A basic helix-loop-helix protein 44) (bHLHa44)
Gene Name	ASCL4 BHLHA44 HASH4
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	121549
Human Swiss-Prot Number	Q6XD76
Alternative Names	
Background	

