



LRIT3 rabbit pAb

Cat No.:ES19046

For research use only

Overview

Product Name	LRIT3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LRIT3 AA range: 156-206
Specificity	This antibody detects endogenous levels of human LRIT3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Leucine-rich repeat, immunoglobulin-like domain and transmembrane domain-containing protein 3
Gene Name	LRIT3
Cellular localization	Cell projection, dendrite . Perikaryon . Endoplasmic reticulum membrane ; Single-pass type I membrane protein . Punctate expression at dendrite tips. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	345193
Human Swiss-Prot Number	Q3SXY7
Alternative Names	
Background	

