



HERC1 rabbit pAb

Cat No.:ES19049

For research use only

Overview

Product Name	HERC1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HERC1 AA range: 509-559
Specificity	This antibody detects endogenous levels of human HERC1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Probable E3 ubiquitin-protein ligase HERC1 (EC 6.3.2.-) (HECT domain and RCC1-like domain-containing protein 1) (p532) (p619)
Gene Name	HERC1
Cellular localization	Membrane; Peripheral membrane protein. Cytoplasm, cytosol. Golgi apparatus. Recruited onto actin-rich surface protrusions.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	8925
Human Swiss-Prot Number	Q15751
Alternative Names	
Background	

