



CHIC1 rabbit pAb

Cat No.:ES19053

For research use only

Overview

Product Name	CHIC1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CHIC1 AA range: 51-101
Specificity	This antibody detects endogenous levels of human CHIC1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Cysteine-rich hydrophobic domain 1 protein (Brain X-linked protein)
Gene Name	CHIC1 BRX
Cellular localization	Cell membrane . Cytoplasmic vesicle . Also present at a Golgi-like vesicular compartment and at scattered vesicles. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	53344
Human Swiss-Prot Number	Q5VXU3
Alternative Names	
Background	

