

WDR81 rabbit pAb

Cat No.:ES19055

For research use only

Overview

Product Name	WDR81 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human WDR81 AA range: 525-575
Specificity	This antibody detects endogenous levels of human WDR81
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	WD repeat-containing protein 81
Gene Name	WDR81
Cellular localization	Early endosome membrane ; Peripheral membrane protein . Late endosome membrane . Lysosome membrane . Cytoplasmic vesicle, autophagosome membrane . Mitochondrion . Cytoplasm, cytosol .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	124997
Human Swiss-Prot Number	Q562E7
Alternative Names	
Background	

