



ENT1 rabbit pAb

Cat No.:ES19057

For research use only

Overview

Product Name	ENT1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ENT1 AA range: 245-295
Specificity	This antibody detects endogenous levels of human ENT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	HERV-T_19q13.11 provirus ancestral Env polyprotein (Envelope polyprotein) (HERV-T Env protein) [Includes: Surface protein (SU); Transmembrane protein (TM)]
Gene Name	
Cellular localization	Cell membrane ; Single-pass membrane protein . At the origin, this retroviral envelope protein was localized in the virion. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	P61550
Alternative Names	
Background	

