

# BPIB6 rabbit pAb

Cat No.:ES19061

For research use only

## Overview

Product Name	BPIB6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human BPIB6 AA range: 268-318
Specificity	This antibody detects endogenous levels of human BPIB6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	BPI fold-containing family B member 6 (Bactericidal/permeability-increasing protein-like 3)
Gene Name	BPIFB6 BPIL3
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	128859
Human Swiss-Prot Number	Q8NFQ5
Alternative Names	
Background	

