

PPM1A rabbit pAb

Cat No.:ES19063

For research use only

Overview

Product Name	PPM1A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PPM1A AA range: 284-334
Specificity	This antibody detects endogenous levels of human PPM1A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein phosphatase 1A (EC 3.1.3.16) (Protein phosphatase 2C isoform alpha) (PP2C-alpha) (Protein phosphatase 1A)
Gene Name	PPM1A PPPM1A
Cellular localization	Nucleus . Cytoplasm, cytosol . Membrane ; Lipid-anchor . Weakly associates at the membrane and N-myristoylation mediates the membrane localization. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	5494
Human Swiss-Prot Number	P35813
Alternative Names	
Background	

