



S35B3 rabbit pAb

Cat No.:ES19066

For research use only

Overview

Product Name	S35B3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human S35B3 AA range: 258-308
Specificity	This antibody detects endogenous levels of human S35B3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Adenosine 3'-phospho 5'-phosphosulfate transporter 2 (3'-phosphoadenosine 5'-phosphosulfate transporter) (PAPS transporter 2) (Solute carrier family 35 member B3)
Gene Name	SLC35B3 C6orf196 PAPST2 CGI-19
Cellular localization	Golgi apparatus membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	51000
Human Swiss-Prot Number	Q9H1N7
Alternative Names	
Background	

