

EBLN2 rabbit pAb

Cat No.:ES19067

For research use only

Overview

Product Name	EBLN2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human EBLN2 AA range: 80-130
Specificity	This antibody detects endogenous levels of human EBLN2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Endogenous Bornavirus-like nucleoprotein 2 (Endogenous Borna-like N element-2) (EBLN-2)
Gene Name	EBLN2 GK006
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	55096
Human Swiss-Prot Number	Q6P2I7
Alternative Names	
Background	

