



MREG rabbit pAb

Cat No.:ES19068

For research use only

Overview

Product Name	MREG rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MREG AA range: 98-148
Specificity	This antibody detects endogenous levels of human MREG
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Melanoregulin (Dilute suppressor protein homolog)
Gene Name	MREG DSU HDCGA21P
Cellular localization	Apical cell membrane ; Peripheral membrane protein . Melanosome membrane ; Lipid-anchor . Lysosome membrane ; Lipid-anchor . Cytoplasmic vesicle membrane . Localizes to the inner segment and basal outer segment of rods in the retina. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	55686
Human Swiss-Prot Number	Q8N565
Alternative Names	
Background	

