



VP33A rabbit pAb

Cat No.:ES19070

For research use only

Overview

Product Name	VP33A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human VP33A AA range: 272-322
Specificity	This antibody detects endogenous levels of human VP33A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Vacuolar protein sorting-associated protein 33A (hVPS33A)
Gene Name	VPS33A
Cellular localization	Cytoplasmic vesicle . Late endosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Lysosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Early endosome . Cytoplasmic vesicle, autophagosome . Cytoplasmic vesicle, clathrin-coated vesicle .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	65082
Human Swiss-Prot Number	Q96AX1
Alternative Names	
Background	

