

AFF2 rabbit pAb

Cat No.:ES19073

For research use only

Overview

Product Name	AFF2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human AFF2 AA range: 1104-1154
Specificity	This antibody detects endogenous levels of human AFF2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	AF4/FMR2 family member 2 (Fragile X E mental retardation syndrome protein) (Fragile X mental retardation 2 protein) (FMR2P) (Protein FMR-2) (Protein Ox19)
Gene Name	AFF2 FMR2 OX19
Cellular localization	Nucleus speckle . When splicing is inhibited, accumulates in enlarged speckles.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	2334
Human Swiss-Prot Number	P51816
Alternative Names	
Background	

