

CLC9A rabbit pAb

Cat No.:ES19081

For research use only

Overview

Product Name	CLC9A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CLC9A AA range: 33-83
Specificity	This antibody detects endogenous levels of human CLC9A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	C-type lectin domain family 9 member A
Gene Name	CLEC9A UNQ9341/PRO34046
Cellular localization	Membrane ; Single-pass type II membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	283420
Human Swiss-Prot Number	Q6UXN8
Alternative Names	
Background	

