

FBX8 rabbit pAb

Cat No.:ES19090

For research use only

Overview

Product Name	FBX8 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FBX8 AA range: 186-236
Specificity	This antibody detects endogenous levels of human FBX8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	F-box only protein 8 (F-box/SEC7 protein FBS)
Gene Name	FBXO8 FBS FBX8 DC10 UNQ1877/PRO4320
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	26269
Human Swiss-Prot Number	Q9NRD0
Alternative Names	
Background	

