



NXN rabbit pAb

Cat No.:ES19092

For research use only

Overview

Product Name	NXN rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NXN AA range: 271-321
Specificity	This antibody detects endogenous levels of human NXN
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Nucleoredoxin (EC 1.8.1.8)
Gene Name	NXN NRX
Cellular localization	Cytoplasm, cytosol . Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	64359
Human Swiss-Prot Number	Q6DKJ4
Alternative Names	
Background	

