



MEX3A rabbit pAb

Cat No.:ES19093

For research use only

Overview

Product Name	MEX3A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MEX3A AA range: 48-98
Specificity	This antibody detects endogenous levels of human MEX3A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	RNA-binding protein MEX3A (RING finger and KH domain-containing protein 4)
Gene Name	MEX3A RKHD4
Cellular localization	Cytoplasm . Nucleus . Cytoplasm, P-body . Predominantly expressed in the cytoplasm and shuttles between the cytoplasm and the nucleus through the CRM1 export pathway.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	92312
Human Swiss-Prot Number	A1L020
Alternative Names	
Background	

