

FRG2B rabbit pAb

Cat No.: ES19095

For research use only

Overview

Product Name FRG2B rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1:1000-2000

Immunogen Synthesized peptide derived from human FRG2B AA

range: 196-246

Specificity This antibody detects endogenous levels of human

FRG2B

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Protein FRG2-like-1 (FSHD region gene 2 protein

family member B) (HSA10-FRG2)

Gene Name FRG2B
Cellular localization Nucleus.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/ml

Observed band

Human Gene ID 441581 Human Swiss-Prot Number Q96QU4

Alternative Names

+86-27-59760950

Background

