

RNF31 rabbit pAb

Cat No.:ES19097

For research use only

Overview

RNF31 rabbit pAb
Rabbit
WB
Human; Mouse
WB 1:1000-2000
Synthesized peptide derived from human RNF31 AA range: 833-883
This antibody detects endogenous levels of human RNF31
Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Store at -20 $^\circ\!\mathrm{C}$. Avoid repeated freeze-thaw cycles.
E3 ubiquitin-protein ligase RNF31 (EC 6.3.2)
(HOIL-1-interacting protein) (HOIP) (RING finger
protein 31) (Zinc in-between-RING-finger
ubiquitin-associated domain protein)
RNF31 ZIBRA
Cytoplasm .
The antibody was affinity-purified from rabbit
antiserum by affinity-chromatography using
epitope-specific immunogen.
Polyclonal
1 mg/ml
55072
Q96EP0



ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C