



FNIP2 rabbit pAb

Cat No.:ES19098

For research use only

Overview

Product Name	FNIP2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FNIP2 AA range: 745-795
Specificity	This antibody detects endogenous levels of human FNIP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Folliculin-interacting protein 2 (FNIP1-like protein) (O6-methylguanine-induced apoptosis 1 protein)
Gene Name	FNIP2 FNIPL KIAA1450 MAPO1
Cellular localization	Lysosome membrane . Cytoplasm . Colocalizes with FLCN in the cytoplasm. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	57600
Human Swiss-Prot Number	Q9P278
Alternative Names	
Background	

