

TBA3E rabbit pAb

Cat No.:ES19100

For research use only

Overview

Product Name	TBA3E rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TBA3E AA range: 320-370
Specificity	This antibody detects endogenous levels of human TBA3E
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Tubulin alpha-3E chain (Alpha-tubulin 3E)
Gene Name	TUBA3E
Cellular localization	Cytoplasm, cytoskeleton .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	112714
Human Swiss-Prot Number	Q6PEY2
Alternative Names	
Background	

