

ZBED1 rabbit pAb

Cat No.:ES19102

For research use only

Overview

Product Name	ZBED1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZBED1 AA range: 329-379
Specificity	This antibody detects endogenous levels of human ZBED1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger BED domain-containing protein 1 (Putative Ac-like transposable element) (dREF homolog)
Gene Name	ZBED1 ALTE DREF KIAA0785 TRAMP
Cellular localization	Nucleus, PML body . Nucleus . In granular structures. . ; Nucleus, PML body. (Microbial infection) Upon interaction with human adenovirus early E1A protein, the protein is redistributed to the peripheral areas of PML bodies. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9189
Human Swiss-Prot Number	O96006
Alternative Names	
Background	

