

## ARPC2 rabbit pAb

Cat No.: ES19103

For research use only

## Overview

Product Name ARPC2 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse; Rat Recommended dilutions WB 1:1000-2000

Immunogen Synthesized peptide derived from human ARPC2 AA

range: 147-197

**Specificity** This antibody detects endogenous levels of human

ARPC2

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Actin-related protein 2/3 complex subunit 2 (Arp2/3

complex 34 kDa subunit) (p34-ARC)

Gene Name ARPC2 ARC34 PRO2446

**Cellular localization** Cytoplasm, cytoskeleton . Cell projection . Cell

junction, synapse, synaptosome. Nucleus.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/ml

**Observed band** 

Human Gene ID 10109 Human Swiss-Prot Number 015144

**Alternative Names** 

+86-27-59760950

**Background** 

