



CAB45 rabbit pAb

Cat No.:ES19104

For research use only

Overview

Product Name	CAB45 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CAB45 AA range: 63-113
Specificity	This antibody detects endogenous levels of human CAB45
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	45 kDa calcium-binding protein (Cab45) (Stromal cell-derived factor 4) (SDF-4)
Gene Name	SDF4 CAB45 PSEC0034
Cellular localization	[Isoform 1]: Golgi apparatus lumen .; [Isoform 5]: Cytoplasm . Cell membrane . Cell projection, bleb. Isoform 5 colocalizes with STX3 and STXBP1 isoform 2 at the plasma membrane and cell surface blebs. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	51150
Human Swiss-Prot Number	Q9BRK5
Alternative Names	
Background	

